# Breast Self-Examination Aid

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## Table of Contents

Liv Breast Self-Examination.....	1
When to Perform a Breast Self-Examination.....	2
Why Self Checking Matters.....	2
Care and Storage.....	2
Instructions.....	3
Important Information.....	5
Specific Use Disclaimer.....	5

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## Liv Breast Self-Examination Aid

**Consistency and Routine** are the two key words in your optimal breast self-exam schedule. The more consistent your approach is to the exam, the easier it will be. Before long, it will become second nature and part of your regular health and wellness routine.

It is important to understand the benefits and limitations of a self-examination. A self-examination is a positive action toward understanding what looks and feels normal for your breasts. Regular breast self-examination does not replace the need for regular doctor visits and screenings, so please make sure to schedule appointments according to your age and risk factor. Use the table below to determine your breast health routine and stick to it.

age	breast health routine
20	Monthly breast self-exams
20-40	Monthly breast self-exams + a yearly clinical breast exam by a health care provider
35-40	A baseline mammogram as a comparison with future mammograms
40+	Monthly breast self-exams + a yearly clinical breast exam by a health care provider + a yearly mammogram

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## When to perform a breast self-examination

### **It is important to examine your breasts around the same day every month!**

If you're menstruating, the best time for your self-exam is two to four days after your period has ended because your breast tissue will be less swollen. If you've reached menopause, or are experiencing irregular periods, pick the same calendar day each month to do your self-examination.

## Why self-checking matters

Throughout a woman's life, her breast tissue goes through many natural changes; such as menstruation, puberty, pregnancy, breast feeding, menopause, and normal weight changes. There can be lumpiness, swelling, and discomfort, unique to each woman. It is important to be aware and know what is normal for you, so that you have a better chance of recognizing changes. Immediately discuss any changes with your doctor. Every woman's breasts are as unique as the woman herself, it is important to know what is normal for oneself.

## Care and Storage

If the Liv Aid becomes soiled, simply use mild soap to wash it, rinse with warm water and let dry. Store the Liv Aid at room temperature. Never put the Aid into the microwave to warm the device. Just hold the device in your hand until it reaches a comfortable body temperature.

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

## Instructions

### First Look

Look into the mirror with your arms at your sides. Do you notice any change in the shape? Any swelling? Do you notice any skin changes such as dimpling or puckering or an "orange peel" look to your skin?


Next, look again for any changes with your arms held over your head. Does it look normal for you?

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

## Then Feel

Now lay down to flatten out your breast tissue. It is best to place a small pillow under your right shoulder. Put your right arm behind your head so that you can easily examine the area under your armpit. Place the Liv Aid over the area of the breast tissue you want to examine first. Glide your fingers over the Liv Aid in a circular motion, pushing down slightly. Once you examine the entire area that is covered by the Liv Aid, move the Aid to the next area. It is important to examine the entire breast area.

Repeat the entire exam 3 times increasing the pressure each time. You should do your exam using light, medium and firm pressure.

Next remove the pillow and put it under your left shoulder. Repeat all of the steps to examine your left breast.

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## Important

It is very important that you understand three necessary components of a comprehensive breast-screening program.

- ♥ Breast self-examination (with or without the Liv Aid)
- ♥ Regular clinical breast exams by your doctor
- ♥ Mammograms

Breast self-examination with or without the Liv Aid is not a replacement for a mammogram or clinical breast exam by your doctor. It is very important that you have regular clinical exams and mammograms as advised by your doctor. Report any changes found either with your bare hand or the Liv Aid exam to your doctor for proper evaluation.

If you discover a lump, puckering or discharge during either the barehanded breast exam or the breast exam using the Liv Aid, it is important that you see your doctor as soon as possible.

**SPECIFIC USE DISCLAIMER:** There is no guarantee that the Liv Aid will detect breast cancer, lumps or any other breast abnormality or disease. Breast self-examination is not a substitute for regular medical examinations and are only a part of good breast health care. Clinical breast examinations and mammograms should be regularly performed by a qualified physician. Changes found during breast self-examination should be promptly reported to your physician for proper evaluation. Women with prior breast health problems or disease should rely primarily on more frequent clinical examinations by a qualified physician. Women with reconstructive surgery of the breast or breast implants should consult with their physician before using this product."

Liv is a trademark of Olivia Newton-John use with permission.

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*English*

## **Remember your breast health routine**

**CROSS YOUR HEART. MAKE A PROMISE.** Remind yourself by writing on the Liv Promise Sticker. Fill in the dates of your monthly breast self-exam, your next doctor appointment and your next mammogram appointment. Peel it off and stick somewhere visible as a reminder, such as on your bathroom mirror, your bed stand, or your bathroom cabinet. The Liv heart will help you to remember your promise to keep a good breast health routine and an active approach in overall health and wellness.

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